DERMA PRENE®
Latex-Free, Powder-Free Surgical Gloves
Gants chirurgicaux sans latex et sans poudre

Ansell

www.ansellcanada.ca

Understanding Latex Allergies

Adverse Skin Reactions Associated with Glove Use

Adverse reaction to NRL (natural rubber latex) gloves can range from mild irritation to a serious allergic response. The four major types of adverse skin reactions associated with NRL glove use are immediate hypersensitivity (Type I or latex allergy) delayed hypersensitivity (Type IV or contact dermatitis), irritant contact dermatitis and glove powder irritation.

Immediate Type I Response: Latex Allergy

A Type I response is a reaction to residual proteins found in latex. While there are more than 250 different types of latex proteins, approximately 20% are allergenic.

The reaction is immediate, typically occurring 5–30 minutes after initial contact. The symptoms commonly include:

- swelling and redness, local to the site of exposure
- non specific symptoms such as itching and burning

The symptoms can spread to areas remote from the site of contact with the glove and be accompanied by the following:

- conjunctivitis
- rhinitis
- bronchial obstruction
- anaphylaxis

Delayed Type IV Response: Allergic Contact Dermatitis

A Type IV allergy is a reaction to specific allergens, such as chemical residues from the glove manufacturing process (commonly chemical accelerators). Three main classes of chemical compounds are used as accelerators: thiurams, dithiocarbamates and mercaptobenzothiazoles (MBTs). Some residual accelerators can be a source of skin irritation. The response is delayed rather than immediate, usually occurring 6–48 hours after initial contact, although symptoms can last for up to four days.

The symptoms include:

- erythema
- swelling
- cracking
- itching
- weeping
- dryness of the skin at the site, although dermatitis may extend beyond the area of contact

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Managing Latex Allergies in Healthcare

Gestion des allergies au latex

The Prevalence of Natural Rubber Latex Allergies

The introduction of universal glove precautions (commonly referred to as “standard precautions”) to prevent the transmission of disease has increased the exposure of healthcare professionals to NRL gloves. This increased level of exposure resulted in an apparent initial rise in the prevalence of NRL allergies, particularly in the 1980s when glove precautions were first implemented. Mandatory glove policies also meant that gloves were in greater demand. This demand often exceeded supply, allowing gloves from substandard manufacturers and highly allergenic gloves to infiltrate the market. The incidence of NRL allergy is now decreasing due to steady improvement in manufacturing technologies, latex allergy education and the development of latex alternatives.

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Staff and Patients

The key to managing latex allergies and adverse glove reactions in healthcare professionals and patients lies in correct recognition and appropriate action. Establishing a latex-free environment, including the use of NRL-free or powder-free policies and procedures, is an important preventative measure to reduce the risk of adverse glove reactions.

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This can be done by forming a multidisciplinary committee responsible for developing uniform policies and procedures to protect both patients and healthcare professionals. The committee should be responsible for developing and maintaining the following:

- a latex-safe standard operating procedure manual
- a latex-safe allergy cart or means for staff to access non-latex items
- a proactive occupational health program
- correct latex product identification
- a policy of removal of latex items from patient care areas

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
Are you protected from BOTH type of allergies?

Today, healthcare professionals are highly conscious of the allergenic risks that are linked to natural or synthetic rubber latex gloves. Among these risks are Type IV allergies, which some clinical evidence shows as representing up to 28% of glove-related allergic reactions, with up to 90% of these Type IV allergies being caused by vulcanization accelerators.*

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Reducing Glove Allergenicity

Ansell is the only major glove manufacturer to offer surgical gloves that are both natural rubber latex-free and accelerator-free

Ansell is committed to reducing glove allergenicity at every stage of the manufacturing process. State-of-the-art manufacturing processes ensure that residual proteins levels are reduced at each step of the production, from the initial preparation of the raw latex to the use of PEARL (Protein and Endogenous Allergen Reduction Leaching) technology when necessary.

Ansell is also working to reduce glove allergenicity by continually refining its use of chemical accelerators. Different accelerator types, concentrations and processes are being used in the manufacture of gloves to develop the best product possible.

Natural Rubber Latex-Free Gloving Materials

Because of allergy concerns, a growing number of healthcare facilities have eliminated or limited latex products, including gloves. The alternative materials for latex surgical gloves include neoprene (polychloroprene) and polyisoprene. The following table compares these natural rubber latex-free materials.

A Comparison of Synthetic Gloving Materials

| Glove Material | Level of Barrier Protection | Elasticity | Puncture Resistance | Fit and Comfort | Chemical Resistance | Economy | Low Allergen Content | Recommended Use |
|--|-----------------------------|------------|---------------------|-----------------|---------------------|---------|----------------------|--|
|  Polyisoprene | Very Good | Excellent | Good | Excellent | Fair | Good | Very Good | Polyisoprene gloves are recommended as an alternative to natural rubber latex because of their similarity to latex with respect to fit, feel and comfort. They are also recommended for individuals with latex allergies. |
|  Neoprene | Excellent | Excellent | Good | Very Good | Excellent | Good | Excellent | Neoprene gloves are recommended as an alternative to natural rubber latex because of their similarity to latex with respect to fit, feel, comfort and barrier protection. They are also recommended for individuals with latex allergies and sensitive skin, as well as Type IV allergies. |

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Safety and Protection

Infection Control

Viral Penetration. Although not required by the FDA, Ansell's premium line of surgical gloves is tested to ensure protection against viral penetration. The Encore® family of gloves meets the viral penetration requirements according to ASTM F1671.

Cuff. Beaded.

Cuff Roll-Down. Ansell's SureFit™ Technology overcomes the "cuff roll-down" typically experienced while wearing impervious surgical gowns designed to prevent fluid penetration. By maintaining the surgical glove securely around the arm, SureFit™ Technology protects the integrity of the sterile field by shielding the gown cuff and ensures the efficacy of the surgical glove as a barrier protection.

Bone Cement. As long as the preparation and use directions on the bone cement are followed, these gloves provide a protective barrier for handling mixed bone-cement product.

Chemotherapy Drugs. These gloves have been tested for use with chemotherapy drugs. Gloves use for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. List of drugs tested against these gloves are available upon request to our Customer Service Representative.

Skin Protection

Does not contain natural rubber latex. Safe for latex-sensitive (Type I) healthcare professional and their patients.

Powder-Free. Help eliminate powder-related complications and sensitivity.

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Comfort and Fit

Donning

Alcohol Resistant. Ansell's proprietary polyurethane inner coating will not break down. It is resistant to alcohol-based hand disinfectants and adds another layer of safety between your skin and the glove.

Damp-Hand Donnability. Enhanced polymer coating allows easy dry and damp-hand donning.

Double Donnability. Easy and comfortable double-gloving.

Hand Fatigue

Shape of Gloves/Former. Produced on new, ergonomically-designed formers to enhance fit and comfort, and to reduce hand fatigue and stress.

Stretchability (Modulus). Wearers experience a more secure fit that provide greater comfort without bagginess.

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In addition to the overall benefits that Ansell Derma Prene® Surgical gloves offer, each glove in this Family provides unique features for specific needs.

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Yip, E.S., Comments to the Maine legislature on proposed prohibition of sale of non-sterile latex gloves. 2003.
O'Gilvie W., Latex sensitization in the healthcare setting. Device Bulletin 9601 ; April 1996, The Medical Devices Agency of the Department of Health.
Susman E., AAAAI: Latex sensitivity infrequent in healthcare workers. In: Doctor's Guide Global Edition. 2003.
Susman G., The effects of interventions and glove changes in health care workers with latex allergy. AnnAllergy Asthma Immunol; 90: 179-80

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® & TM are trademarks owned by Ansell Healthcare Products LLC or one of its affiliates. Printed in Canada. ©2009 All Rights Reserved.