

Encore®

Encore® surgical gloves are powder-free and made from natural rubber latex. The gloves have a textured finish for superior grip and handling, and their low allergen levels help reduce the risk of latex sensitization in healthcare workers and patients.

Safety & Protection

Infection Control

Viral Penetration. Although not an FDA requirement, Ansell's premium line of surgical gloves is tested to ensure protection against viral permeation. Encore® has been tested and passed the viral penetration according to ASTM F1671-97a.

Cuff Roll-Down. This glove has Ansell's patent-pending SureFit Technology™, which provides surgeons, nurses and other health-care professionals with superior barrier protection during operating room procedures.

Bone Cement. As long as the preparation and use directions on the bone cement are followed, this glove provides a protective barrier for handling mixed bone-cement product.

Embryo Use. This glove is appropriate and accepted for embryo use.

Skin Protection

Low Protein Content. This glove contains 50 µg/g or less of water/extractable protein per gram.

Caution: Safe use of these gloves by latex-sensitized individuals has not been established.

Comfort & Fit

Donning

Alcohol Resistant. Ansell's proprietary polyurethane inner coating will not break down. It is resistant to alcohol-based hand disinfectants and adds another layer of safety between your skin and the latex.

Damp-Hand Donnability. Polymer coating for dry- and damp-hand donnability.

Double Donnability. This glove is recommended to be an over glove because of its textured surface.

Hand Fatigue

Shape of Glove/Former. Designed with Anatomic Curved Former for a more natural fit and less fatigue. This exclusive Ansell former is the result of a staff dedicated to the development and continuous improvement of formers.

Stretchability (Modulus). Because this glove is latex, wearers experience a more secure fit that provides greater comfort without being baggy.

Tactile Sensitivity

Dependability. The thickness of these gloves is optimized for the best balance of protection and tactile sensitivity, all while allowing for a comfortable fit during the majority of surgical procedures.

Grip

Application Driven. A textured finish means unequalled reliability when using instruments that require greater gripping capability.



ENCORE®



Ansell

Everything you touch...we touch.