

# Encore<sup>®</sup> MicrOptic<sup>®</sup>

Encore<sup>®</sup> MicrOptic<sup>®</sup> latex, powder-free surgical gloves offer exceptional dexterity, tactile sensitivity, and at the same time, outstanding strength and barrier protection. The Encore<sup>®</sup> MicrOptic<sup>®</sup> is ideal for delicate procedures in ophthalmic, micro, and cardiothoracic surgery. Its micro-thickness allows for comfortable double-gloving, and its brown color reduces glare from surgical lights.

## Safety & Protection

### Infection Control

**Viral Penetration.** Although not an FDA requirement, Ansell's premium line of surgical gloves is tested to ensure protection against viral permeation. Encore<sup>®</sup> has been tested and passed the viral penetration according to ASTM F1671-97a.

**Cuff Roll-Down.** This glove has Ansell's patent-pending SureFit Technology<sup>™</sup>, which provides surgeons, nurses and other health-care professionals with superior barrier protection during operating room procedures.

**Bone Cement.** As long as the preparation and use directions on the bone cement are followed, this glove provides a protective barrier for handling mixed bone-cement product.

**Embryo Use.** This glove is appropriate and accepted for embryo use.

### Skin Protection

**Low Protein Content.** This glove contains 50 µg/g or less of water/extractable protein per gram.

Caution: Safe use of these gloves by latex-sensitized individuals has not been established.

## Comfort & Fit

### Donning

**Alcohol Resistant.** Ansell's proprietary polyurethane inner coating will not break down. It is resistant to alcohol-based hand disinfectants and adds another layer of safety between your skin and the latex.

**Damp-Hand Donnability.** Polymer coating for dry- and damp-hand donnability.

**Double Donnability.** This glove can be used as an under glove or an over glove depending upon the user's preference.

### Hand Fatigue

**Shape of Glove/Former.** Designed with Anatomic Curved Former for a more natural fit and less fatigue. This exclusive Ansell former is the result of a staff dedicated to the development and continuous improvement of formers.

**Stretchability (Modulus).** Because this glove is latex, wearers experience a more secure fit that provides greater comfort without being baggy.

## Tactile Sensitivity

**Dependability.** This glove is recommended for surgical procedures requiring a high level of tactile sensitivity.

## Grip

**Application Driven.** Anti-slip smooth finish prevents slip.



MICROPTIC<sup>®</sup>

Ansell

Everything you touch...we touch.