

Material Natural rubber from latex in compliance with ASTM D 1076

Packaging 1 Sterile Pair/Polyethylene Package – 50 Pairs/Dispenser Box – 4 Boxes/Case

Product Numbers

| SIZE (identified inside glove) | ASPN (Ansell Standard Product Number) |
|-----------------------------------|--|
| 5 1/2 | 5787000 |
| 6 | 5787001 |
| 6 1/2 | 5787002 |
| 7 | 5787003 |
| 7 1/2 | 5787004 |
| 8 | 5787005 |
| 8 1/2 | 5787006 |
| 9 | 5787007 |

Dimensional Specifications (minimum length 295mm)

| SIZE | Minimum Length mm (inches) | Average Palm Width mm (inches), ±4mm |
|-------|-------------------------------|---|
| 5 1/2 | 296 (11.7) | 70 (2.7) |
| 6 | 297 (11.7) | 76 (3.0) |
| 6 1/2 | 298 (11.7) | 82 (3.2) |
| 7 | 298 (11.7) | 89 (3.5) |
| 7 1/2 | 298 (11.7) | 95 (3.7) |
| 8 | 300 (11.8) | 102 (4.0) |
| 8 1/2 | 300 (11.8) | 108 (4.3) |
| 9 | 303 (11.9) | 115 (4.5) |

Thickness (single)

| Data from sampling | Average Thickness | Range |
|----------------------|-------------------|--------------|
| Fingertip, mm (mils) | 0.18 (7.1) | ±0.018 (0.7) |
| Cuff, mm (mils) | 0.11 (5) | ±0.011 (0.5) |

Physical Properties (ASTM, EN and ISO Requirements)

| | Before Aging | After Aging |
|---|--------------|--------------|
| Freedom from Holes (inspection level I) | 1 AQL* | N/A |
| Tensile, MPa (p.s.i.) minimum | 24.0 (3,480) | 18.0 (2,610) |
| Elongation (%) minimum | 750 | 560 |
| Modulus, MPa (p.s.i.) maximum | 5.5 (800) | N/A |

Manufacturing specifications meet or exceed all the above.

| | |
|--------------------------|--|
| Color | Brown |
| Accelerators | Zinc dibutyldithiocarbamate (ZDBC) Zinc diethyldithiocarbamate (ZDEC) |
| Powder-Free | Eliminates powder-related post-operative complications and sensitivity |
| Certification | Manufacturing Facility Registered to ISO 9001 and 13485 (The latter supersedes EN 46001) Contains 50 micrograms or less of total water-extractable protein per gram |
| Quality Control | In compliance with ASTM D 3577-00, Type I |
| Sterilization | Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000 |
| References | ANSI/ASQC Z1.4 Sampling Procedures and Tables for Inspection by Attributes |
| Country of Origin | Sri Lanka |
| Storage | Cool environment, shielded from ozone and UV light |

*Acceptable Quality Level

Except for applicable ASTM requirements, the above information represents typical data for the product and should not be construed or used as an actual product performance specification.

For information and samples, call Ansell Healthcare at 800-952-9916, or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
Ansell, Encore and MicroOptic are registered trademarks owned by Ansell Healthcare Products LLC or one of its affiliates. ©2009 All Rights Reserved.

SS240 04/2009



Ansell
Everything you touch...we touch.