

MICRO-TOUCH® Nitra-Tex® Sterile

Micro-Touch® Nitra-Tex® Sterile exam gloves feature a 12" extended cuff and textured fingertips, providing superior barrier protection and grip in a sterile, latex-free exam glove.

Safety & Protection

Infection Control

Allergy Concerns. This examination glove is made of 100% nitrile. It does not contain natural rubber latex thus addressing Type I (latex) allergy concerns.

Sterilization. Gamma sterilized minimum 25kGy to provide sterility.

Quality Benchmarks. Meets or exceeds ASTM examination glove standards. Manufactured within the quality guidelines of ISO 13485:2003, ISO 9001:2008 and FDA-QSR. FDA approved for handling chemotherapy drugs.¹

Cuff. Extended 12" beaded cuff for extra protection.

Resistance. Nitrile film has excellent puncture and chemical resistance. Nitrile gloves are recommended for conditions where high strength and chemical protection are required.

This glove is tested for use with the following 8 chemotherapy drugs: Vincristine Sulfate, Cyclophosphamide, Doxorubicin Hydrochloride, 5-Fluorouracil, Methotrexate, Etoposide, Paclitaxel and ThioTEPA.

Skin Protection

Powder-Free. To help eliminate powder-related complications and sensitivity.

Comfort & Fit

Hand Fatigue

Stretchability (Modulus). Because this glove is nitrile wearers experience a good fit and feel due to high elasticity and memory.

Shape of Glove/Former. Ambidextrous.

Tactile Sensitivity

Dependability. This glove is recommended for sterile exam (non-surgical) procedures, where a latex-free exam glove is required or preferred.

Memory. Nitrile has high elasticity and good memory, allowing the film to adapt to the wearer's hand.

Grip

Application Driven. A textured fingertip finish provides extra protection and reliability during sterile exam (non-surgical) procedures.

Ordering Information

SIZE	Pairs Product#	Singles Product#
Small	6034151	6034251
Medium	6034152	6034252
Large	6034153	6034253
X-Large	6034154	6034254

Pairs packed 50 pairs per dispenser box – 4 boxes per case

Singles packed 100 gloves per dispenser box – 4 boxes per case

1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.



MICRO-TOUCH®
NITRA-TEX® STERILE

Ansell

Everything you touch...we touch.