

Material 100% polychloroprene (neoprene); no natural rubber

Dimensional Specifications

SIZE	Minimum Length mm (inches)	Average Palm Width mm (inches)
X-Small	240 (9.4)	76 (3.0)
Small	240 (9.4)	85 (3.3)
Medium	240 (9.4)	95 (3.7)
Large	240 (9.4)	106 (4.2)
X-Large	240 (9.4)	115 (4.5)

Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.16 (6.3)
Palm, mm (mils)	0.13 (5.1)
Cuff, mm (mils)	0.11 (4.3)

Physical Properties (ASTM Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level 1)	2.5 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	14.0 (2,031)	14.0 (2,031)
Elongation (%) minimum	700	600

Ansell manufacturing specifications meet or exceed ASTM specifications.

Color	Green
Accelerators	Zinc dibutylthiocarbamate (ZDBC) – Residual Concentration: Not Detectable
Certification	Manufacturing facility registered to ISO 9001 FDA approved for handling chemotherapy drugs ¹
Donning Agent	Polymer coating
Residual Powder Content	<2.0 mg/glove
Quality Control	In compliance with ASTM D6977 examination gloves standards
References	ANSI/ASQ Z1.4 sampling procedures and tables for inspection by attributes
Country of Origin	Indonesia
Storage	Cool environment, shielded from ozone and UV light

*Acceptable quality level – maximum % of gloves with defects

¹ Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.

For additional information and samples, call Ansell at 800-952-9916 or visit us at www.ansellhealthcare.com


Ansell Healthcare Products LLC | ansellhealthcare.com

200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)

Ansell, Micro-Touch and Affinity are registered trademarks owned by Ansell Limited or one of its affiliates. ©2011 Ansell Healthcare Products LLC. All Rights Reserved.

SS270 6/2011

Printed on recycled paper using soy inks.



Ansell
Everything you touch...we touch.