

# MICRO-TOUCH® Elite®

MICRO-TOUCH®  
ELITE®

Micro-Touch® Elite® powder-free, stretch vinyl examination gloves are produced with a polyvinyl chloride formulation that provides a more comfortable fit than traditional vinyl gloves. They are an excellent alternative for latex-sensitive individuals for use in low-risk clinical applications that do not involve blood-borne pathogens, viruses or bodily fluids. This glove is not formulated with DEHP.

## Safety & Protection

### Infection Control

**Allergy Concerns.** This examination glove is made from polyvinyl chloride. It does not contain natural rubber latex, thus addressing Type I (latex) allergy concerns. It also does not contain chemical curing agents.

**Sterilization.** This is a non-sterile examination glove.

**Quality Benchmarks.** Meets or exceeds ASTM examination glove standards. Manufactured within the quality guidelines of ISO 13485:2003, ISO 9001:2008 and FDA-QSR.

**Cuff.** Inverted, 9.5" beaded cuff.

### Skin Protection

**Powder-Free.** To help eliminate powder-related complications and sensitivity.



## Comfort & Fit

### Hand Fatigue

**Stretchability (Modulus).** Stretch vinyl examination gloves are produced with a polyvinyl chloride formulation that provides a more comfortable and secure fit than traditional vinyl gloves.

**Shape of Glove/Former.** Ambidextrous.

## Grip

**Application Driven.** Smooth finish for better tactile sensitivity.

## Ordering Information

SIZE	Product Number
Small	3091
Medium	3092
Large	3093
X-Large	3094

100 Gloves/Dispenser Box – 10 Boxes/Case



**Ansell**

Everything you touch...we touch.