

Material 100% polyvinyl chloride (PVC), no natural rubber

Dimensional Specifications

SIZE	Minimum Length mm (inches)	Average Palm Width mm (inches)
Small	230 (9.1)	85 (3.3)
Medium	230 (9.1)	95 (3.7)
Large	230 (9.1)	105 (4.1)
X-Large	240 (9.4)	115 (4.5)

Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.10 (3.9)
Palm, mm (mils)	0.11 (4.3)

Physical Properties (ASTM Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level 1)	2.5 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	11.0 (1,595)	11.0 (1,595)
Elongation (%) minimum	350	300

Ansell manufacturing specifications meet or exceed ASTM specifications.

Color	Cream
Accelerators	None used in the manufacturing process
Certification	Manufacturing facility registered to ISO 9001 and ISO 13485
Donning Agent	None (powder-free)
Residual Powder Content	<2.0 mg/glove
Quality Control	In compliance with ASTM D5250
References	ANSI/ASQ Z1.4 Sampling Procedures and Tables for Inspection by Attributes
Country of Origin	China
Storage	Cool environment, shielded from ozone and UV light

*Acceptable quality level – maximum % of gloves with defects

**For additional information and samples, call Ansell
at 800-952-9916 or visit us at www.ansellhealthcare.com**

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
Ansell, Micro-Touch and Elite are registered trademarks owned by Ansell Limited or one of its affiliates. ©2011 Ansell Healthcare Products LLC. All Rights Reserved.

SS280 11/2011

Printed on recycled paper using soy inks.



Ansell
Everything you touch...we touch.