

**Material** 100% nitrile rubber (acrylonitrile / butadiene / methacrylic acid terpolymer); no natural rubber or sulfur-based chemical accelerators\*

### Dimensional Specifications

SIZE	Minimum Length mm (inches)	Average Palm Width mm (inches)
X-Small	240 (9.4)	75 (3.0)
Small	240 (9.4)	85 (3.3)
Medium	240 (9.4)	95 (3.7)
Large	240 (9.4)	106 (4.2)
X-Large	240 (9.4)	116 (4.6)

### Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.08 (3.1)
Palm, mm (mils)	0.06 (2.4)
Cuff, mm (mils)	0.05 (2.0)

### Physical Properties (ASTM Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level 1)	2.5 AQL**	N/A
Tensile, MPa (p.s.i.) minimum	14.0 (2,031)	14.0 (2,031)
Elongation (%) minimum	500	400

Ansell manufacturing specifications meet or exceed ASTM specifications.

<b>Color</b>	Pink
<b>Accelerators</b>	None
<b>Certification</b>	Manufacturing facility registered to ISO 9001 and ISO 13485 FDA approved for handling chemotherapy drugs <sup>1</sup>
<b>Donning Agent</b>	None (powder-free)
<b>Residual Powder Content</b>	<2.0 mg/glove
<b>Quality Control</b>	In compliance with ASTM D6319 examination glove standards
<b>References</b>	ANSI/ASQ Z1.4 Sampling Procedures and Tables for Inspection by Attributes
<b>Country of Origin</b>	Malaysia
<b>Storage</b>	Cool environment, shielded from ozone and UV light

\*Does not contain any sulfur-based accelerators such as thiurams, carbamates, thiazoles, guanidine or thioureas

\*\*Acceptable quality level – maximum % of gloves with defects

1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.

**For additional information and samples, call Ansell at 800-952-9916 or visit us at [www.ansellhealthcare.com](http://www.ansellhealthcare.com)**

Ansell Healthcare Products LLC | [ansellhealthcare.com](http://ansellhealthcare.com)  
 200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)  
 Ansell, Micro-Touch and NitraFree are trademarks owned by Ansell Limited or one of its affiliates. ©2011 Ansell Healthcare Products LLC. All Rights Reserved.  
 The COLOR PINK® and PINK NITRILE® are registered trademarks owned by Colur World, LLC. 2006 All Rights Reserved.

SS460 6/2011

Printed on recycled paper using soy inks.



**Ansell**

Everything you touch...we touch.