

**Material** Natural rubber latex, in compliance with ASTM D1076

### Dimensional Specifications

SIZE	Minimum Length mm (inches)	Average Palm Width mm (inches)
<b>X-Small</b>	241 (9.5)	75 (3.0)
<b>Small</b>	241 (9.5)	89 (3.5)
<b>Medium</b>	241 (9.5)	96 (3.8)
<b>Large</b>	241 (9.5)	106 (4.2)
<b>X-Large</b>	241 (9.5)	119 (4.7)

### Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.18 (7.1)
Palm, mm (mils)	0.17 (6.7)

### Physical Properties (ASTM Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level 1)	2.5 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	18.0 (2,611)	14.0 (2,031)
Elongation (%) minimum	650	500

*Ansell manufacturing specifications meet or exceed ASTM specifications.*

<b>Color</b>	Light green
<b>Accelerators</b>	Zinc dibutyldithiocarbamate (ZDBC) Zinc diethyldithiocarbamate (ZEBC)
<b>Certification</b>	Manufacturing facility registered to ISO 9001, EN 46001 and ISO 13485
<b>Donning Agent</b>	None (powder-free)
<b>Residual Powder Content</b>	<2.0 mg/glove
<b>Quality Control</b>	In compliance with ASTM D3578 examination glove standards
<b>References</b>	ANSI/ASQ Z1.4 Sampling Procedures and Tables for Inspection by Attributes
<b>Country of Origin</b>	Thailand
<b>Storage</b>	Cool environment, shielded from ozone and UV light

\*Acceptable quality level – maximum % of gloves with defects

**For additional information and samples, call Ansell at 800-952-9916 or visit us at [www.ansellhealthcare.com](http://www.ansellhealthcare.com)**

Ansell Healthcare Products LLC | [ansellhealthcare.com](http://ansellhealthcare.com)

200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)

Ansell, Micro-Touch and NextStep are registered trademarks owned by Ansell Limited or one of its affiliates. ©2011 Ansell Healthcare Products LLC. All Rights Reserved.