

MICRO-TOUCH® NextStep®

Micro-Touch® NextStep® gloves with a glycerol moisturizer provide superior protection and tactile sensitivity, along with the comfort and advantages of latex. Designed with the healthcare professional's daily work routine in mind, the light green exam gloves are easy to don, textured for reliable grip and nourish skin during use.

Safety & Protection

Infection Control

Allergy Concerns. Advanced manufacturing processes reduce total water-extractable protein to 50 µg or less per gram. Safe use of these gloves by or on latex-sensitized individuals has not been established.

Sterilization. This is a non-sterile examination glove.

Quality Benchmarks. Meets or exceeds ASTM examination glove standards. Manufactured within the quality guidelines of ISO 13485:2003, FDA-QSR and CE Certification.

Cuff. Inverted, 9.5" beaded cuff.

Resistance. Latex remains the gold standard for hand barrier protection due to its strength, proven barrier protection, elasticity, fit, feel and comfort.



Skin Protection

Powder-Free. To help eliminate powder-related complications and sensitivity. This glove also contains a glycerol coating that is designed to hydrate and moisturize skin that has had continual glove-wearing and frequent contact with astringent hand cleansers and disinfectants.

Comfort & Fit

Hand Fatigue

Stretchability (Modulus). Latex has excellent elasticity allowing the glove to stretch easily and return to its original shape.

Shape of Glove/Former. Ambidextrous.

Tactile Sensitivity

Dependability. Because this glove is latex, wearers experience a more secure fit that provides greater comfort without bagginess.

Memory. Latex has excellent elasticity and memory.

Grip

Application Driven. The fully textured finish of the glove provides extra protection for wet or dry grip and against slippage.

Ordering Information

SIZE	Product Number
X-Small	3200
Small	3201
Medium	3202
Large	3203
X-Large*	3204

100 Gloves/Dispenser Box – 10 Boxes/Case

*80 Gloves/Dispenser Box – 10 Boxes/Case



Caution: This product contains natural rubber latex which may cause allergic reactions.

MICRO-TOUCH®
NEXTSTEP®



Ansell

Everything you touch...we touch.