

Material 100% nitrile rubber (acrylonitrile/butadiene/methacrylic acid terpolymer); no natural rubber

Dimensional Specifications

SIZE	Minimum Length mm (inches)	Average Palm Width mm (inches)
X-Small	240 (9.40)	75 (3.0)
Small	240 (9.40)	85 (3.3)
Medium	240 (9.40)	95 (3.7)
Large	240 (9.40)	106 (4.2)
X-Large	240 (9.40)	116 (4.6)

Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.10 (3.9)
Palm, mm (mils)	0.09 (3.5)
Cuff, mm (mils)	0.08 (3.1)

Physical Properties (ASTM, EN, and ISO Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level 1)	2.5 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	14.0 (2,030)	14.0 (2,030)
Elongation (%) minimum	500	400

Ansell manufacturing specifications meet or exceed ASTM specifications.

Color	Blue
Accelerators	Zinc dibutyldithiocarbamate (ZDBC) Zinc diethyldithiocarbamate (ZDEC)
Certification	Manufacturing facility registered to ISO 9001 and ISO 13485 FDA approved for handling chemotherapy drugs ¹
Donning Agent	None (powder-free)
Residual Powder Content	<2.0 mg/glove
Quality Control	In compliance with ASTM D6319 examination glove standards
References	ISO 2859 Sampling Procedures and Tables for Inspection by Attributes
Country of Origin	Malaysia
Storage	Cool environment, shielded from ozone and UV light

Except for applicable ASTM requirements, the above information represents typical data for the product and should not be construed or used as an actual product performance specification.

*Acceptable quality level – maximum % of gloves with defects

1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.

For additional information and samples, call Ansell Healthcare at 800-952-9916 or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
Ansell and Micro-Touch are registered trademarks owned by Ansell Limited or one of its affiliates. ©2010 All Rights Reserved.