

Material 100% nitrile rubber (acrylonitrile/butadiene/methacrylic acid terpolymer); no natural rubber

Dimensional Specifications

SIZE	Minimum Length mm (inches)	Average Palm Width mm (inches)
X-Small	295 (11.6)	75 (3.0)
Small	295 (11.6)	85 (3.3)
Medium	295 (11.6)	95 (3.7)
Large	295 (11.6)	106 (4.2)
X-Large	295 (11.6)	116 (4.6)

Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.12 (4.7)
Palm, mm (mils)	0.10 (3.9)
Cuff, mm (mils)	0.08 (3.1)

Physical Properties (ASTM Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level 1)	2.5 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	14.0 (2,031)	14.0 (2,031)
Elongation (%) minimum	500	400

Ansell manufacturing specifications meet or exceed ASTM specifications.

Color	Blue
Accelerator	Zinc dibutyldithiocarbamate (ZDBC) Zinc diethyldithiocarbamate (ZDEC)
Certification	Manufacturing facility registered to ISO 9001 and ISO 13485 FDA approved for handling chemotherapy drugs ¹
Donning Agent	None (powder-free)
Residual Powder Content	<2.0 mg/glove
Quality Control	In compliance with ASTM D6319 examination glove standards
References	ANSI/ASQ Z1.4 Sampling Procedures and Tables for Inspection by Attributes
Country of Origin	Malaysia
Storage	Cool environment, shielded from ozone and UV light

*Acceptable quality level – maximum % of gloves with defects


1. Specific chemical resistance data available upon request: 1-800-952-9916 (U.S.A. only). Gloves used for protection against chemotherapy drugs must be selected specifically for the type of chemicals used. Review material safety data sheets for the chemicals being used to determine the required level of protection.

For additional information and samples, call Ansell at 800-952-9916 or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
Ansell, Micro-Touch and E.P. are registered trademarks owned by Ansell Limited or one of its affiliates. ©2011 Ansell Healthcare Products LLC. All Rights Reserved.

SS400 4/2011

Printed on recycled paper using soy inks.



Ansell
Everything you touch...we touch.