

# The Original Perry® Style 42®

The Original Perry® Style 42® powdered surgical gloves provide uncompromising quality and durability. Its satin-smooth finish provides maximum tactile sensitivity and is produced with advanced technology to lower powder and protein levels.

## Safety & Protection

### Infection Control

**Viral Penetration.** Although not an FDA requirement, Ansell's line of surgical gloves is tested to ensure protection against viral permeation. Perry 42® has been tested and passed the viral penetration according to ASTM F1671-97a.

**Cuff Roll-Down.** Inverted, beaded cuff to help prevent cuff roll-down.

**Bone Cement.** Always follow manufacturers' instructions when mixing and handling bone cement. This glove provides a protective barrier for handling and mixing bone-cement product. Double gloving is the usual practice.

### Skin Protection

**Low Protein Content.** This glove contains 60 µg/g or less of water/extractable protein per gram. Caution: Safe use of these gloves by latex-sensitized individuals has not been established.



## Comfort & Fit

### Donning

**Alcohol Resistant.** Ansell latex gloves provide adequate barrier protection for splash protection and/or short duration of alcohol contact.

**Damp-Hand Donnability.** USP absorbable dusting powder aids in smooth donning.

**Double Donnability.** USP absorbable dusting powder also aids in smooth double donning.

### Hand Fatigue

**Shape of Glove/Former.** Produced on curved finger formers, ergonomically designed to enhance fit and comfort, and to reduce hand fatigue and stress.

**Stretchability (Modulus).** Latex has excellent elasticity allowing the glove to stretch easily and return to its original shape. Unique to latex is the ability to reseal small holes if breached.

## Tactile Sensitivity

**Dependability.** This glove is recommended for all levels of risk (high, medium, and low) where a sterile, fitted surgeon's glove would be applicable.

## Grip

**Application Driven.** Anti-slip, smooth finish for a comfortable, secure grip.

*Caution: This product contains natural rubber latex, which may cause allergic reactions.*

THE ORIGINAL  
PERRY® STYLE 42®



**Ansell**

Everything you touch...we touch.