

Technical Data

The Original Perry® Style 42® Latex, Powdered Surgical Gloves, 571110x

Material Natural rubber from latex in compliance with ASTM D 1076

Packaging 1 Sterile Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case (200 Pairs/Case)

Product Numbers

SIZE (identified inside glove)	ASPN (Ansell Standard Product Number)
5 1/2	5711100
6	5711101
6 1/2	5711102
7	5711103
7 1/2	5711104
8	5711105
8 1/2	5711106
9	5711107

Dimensional Specifications (minimum length 276mm)

SIZE	Average Length mm (inches)	Average Palm Width mm (inches)
5 1/2	280 (11.0)	70 (2.8)
6	285 (11.2)	76 (3.0)
6 1/2	285 (11.2)	82 (3.2)
7	298 (11.7)	89 (3.5)
7 1/2	298 (11.7)	95 (3.7)
8	298 (11.7)	102 (4.0)
8 1/2	305 (12.0)	109 (4.3)
9	305 (12.0)	115 (4.5)

Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.23 (9.1)
Palm, mm (mils)	0.21 (8.3)
Cuff, mm (mils)	0.17 (6.5)

Physical Properties (ASTM, EN, and ISO Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level I)	1 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	24 (3,480)	18 (2,610)
Elongation (%) minimum	750	560
Modulus, MPa (p.s.i.) maximum	5.5 (800)	N/A

Ansell manufacturing specifications meet or exceed all the above.

Color	White
Accelerators	Zinc diethyldithiocarbamate (ZDEC)
Donning Agent	Absorbable dusting powder, USP
Residual Powder Content	<15 mg/dm ²
Certification	Manufacturing facility registered to ISO 9001, ISO 13485, and EN 46001 Contains 60 micrograms or less of total water-extractable protein per gram per glove
Quality Control	In compliance with ASTM D 3577, Type 1
Sterilization	Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000
References	ANSI/ASQC Z1.4 Sampling Procedures and Tables for Inspection by Attributes
Country of Origin	Sri Lanka
Storage	Cool, dry environment, shielded from ozone and UV light

*Acceptable Quality Level

Except for applicable ASTM requirements, the above information represents typical data for the product and should not be construed or used as an actual product performance specification.



For additional information and samples, call Ansell Healthcare at 800-952-9916 or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
Ansell, Perry and Style 42 are trademarks owned by Ansell Healthcare Products LLC or one of its affiliates. ©2009 All Rights Reserved.

SS350 04/2009

Ansell
Everything you touch...we touch.