

MICRO-TOUCH® Powder



MICRO-TOUCH®
POWDER

Micro-Touch® Powder gloves feature an anti-slip smooth finish for better tactile sensitivity and a straight cuff. This glove is made from natural rubber latex for superior hand barrier protection and is powdered with absorbable dusting powder. Now featured in a 150 count box to help you go GREEN. With more gloves in a box it means less packaging waste.

Safety & Protection

Infection Control


Allergy Concerns. This examination glove contains natural rubber latex which may cause allergic reactions.

Sterilization. This is a non-sterile examination glove.

Quality Benchmarks. Meets or exceeds ASTM examination glove standards. Manufactured within the quality guidelines of ISO 13485:2003, ISO 9001:2008, ISO 4001 and FDA-QSR.

Cuff. Straight, 9.5" cuff.

Resistance. Latex remains the gold standard for hand barrier protection due to its strength, proven barrier protection, elasticity, fit, feel and comfort.



Comfort & Fit

Hand Fatigue

Stretchability (Modulus). Natural rubber latex has excellent elasticity allowing the glove to stretch easily and return to its original shape.

Shape of Glove/Former. Ambidextrous. This glove is manufactured utilizing Ansell's unique ergonomic and concave formers which provide unparalleled fit.

Tactile Sensitivity

Dependability. Because this glove is made of natural rubber latex, wearers experience an excellent fit that provides greater comfort without bagginess.

Memory. Latex has excellent elasticity and memory.

Grip

Application Driven. A smooth finish allows for better tactile sensitivity when performing non-sterile examination procedures.

Ordering Information

SIZE	Product Number
Small	6116301
Medium	6116302
Large	6116303

150 Gloves/Dispenser Box – 10 Boxes/Case



Caution: This product contains natural rubber latex, which may cause allergic reactions.

Ansell

Everything you touch...we touch.