

Derma Prene® IsoTouch™

The Derma Prene® IsoTouch™ surgical glove is powder-free and made from synthetic polyisoprene. The glove delivers all the performance characteristics of natural rubber latex without the risk of natural rubber latex sensitization. The Derma Prene® IsoTouch™ is light green in color for instant identification as a natural rubber latex-free glove, which helps to reduce confusion and risk to both healthcare workers and patients with sensitive skin.

Safety & Protection

Infection Control

Viral Penetration. Although not required by the FDA, Ansell's premium line of surgical gloves is tested to ensure protection against viral permeation. Derma Prene® IsoTouch™ has been tested and passed the viral penetration according to ASTM F1671 and F1670.

Cuff Roll-Down. This glove features Ansell's patent-pending SureFit Technology™. SureFit Technology™ provides surgeons, nurses and other healthcare professionals with superior protection during operating room procedures.

Bone Cement. As long as the preparation and use directions on the bone cement are followed, this glove provides a protective barrier for handling mixed bone-cement product.

Skin Protection

Does not contain natural rubber latex. Made of 100% synthetic polyisoprene and is safe for latex-sensitive (Type I) healthcare professionals and patients.

Comfort & Fit

Donning

Alcohol Resistant. Ansell's proprietary polyurethane inner coating will not break down. It is resistant to alcohol-based hand disinfectants.

Damp-Hand Donnability. Enhanced polymer coating allows easy dry- and damp-hand donning.

Double Donnability. Design and thickness allow for easy and comfortable double-gloving.

Hand Fatigue

Shape of Glove/Former. Produced on new, ergonomically-designed formers to enhance fit and comfort, and to reduce hand fatigue and stress.

Stretchability (Modulus). Because this glove is made from non-latex polyisoprene, wearers experience a more secure fit that provides greater comfort without bagginess.

Tactile Sensitivity

Dependability. The synthetic polyisoprene material delivers all the performance characteristics of natural rubber latex without the risk of latex sensitization.

Grip

Application Driven. Anti-slip smooth finish prevents slippage.



DERMA PRENE®
ISOTOUCH™

Ansell

Everything you touch...we touch.