

**Material** 100% synthetic polyisoprene; no natural rubber

**Packaging** 1 Sterile Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case

### Product Numbers

SIZE (identified inside glove)	ASPN (Ansell Standard Product Number)
5 1/2	20685255
6	20685260
6 1/2	20685265
7	20685270
7 1/2	20685275
8	20685280
8 1/2	20685285
9	20685290

### Dimensional Specifications (minimum length 295mm)

SIZE	Average Length mm (inches)	Average Palm Width mm (inches)
5 1/2	306 (12.0)	74 (2.9)
6	308 (12.1)	78 (3.1)
6 1/2	308 (12.1)	88 (3.5)
7	312 (12.3)	91 (3.6)
7 1/2	311 (12.2)	97 (3.8)
8	311 (12.2)	103 (4.1)
8 1/2	312 (12.3)	109 (4.3)
9	312 (12.3)	115 (4.5)

### Thickness (single)

Data from sampling	Average Thickness
Fingertip, mm (mils)	0.23 (9.1)
Palm, mm (mils)	0.22 (8.7)
Cuff, mm (mils)	0.18 (7.1)

### Physical Properties (ASTM, EN and ISO Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level I)	1 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	17 (2,465)	12 (1,740)
Elongation (%) minimum	650	490
Modulus, MPa (p.s.i.) maximum	7 (1,015)	N/A

Manufacturing specifications meet or exceed all the above.

<b>Color</b>	Green
<b>Accelerators</b>	Zinc/Sodium Dithiocarbamate Zinc Dibutyldithiocarbamate (ZDBC) Zinc 2-Mercaptobenzothiazole
<b>Powder Free</b>	Eliminates powder-related post-operative complications and sensitivity <2mg/glove (Typical Average Value, 0.28mg/glove)
<b>Certification</b>	Manufacturing Facility Registered to ISO 9001, ISO 13485/EN552
<b>Quality Control</b>	In compliance with ASTM D 3577, Type II
<b>Sterilization</b>	Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000
<b>References</b>	ANSI/ASQ Z1.4 Sampling Procedures and Tables for Inspection by Attributes
<b>Country of Origin</b>	Sri Lanka
<b>Storage</b>	Cool environment, shielded from ozone and UV light

\*Acceptable Quality Level

Except for applicable ASTM requirements, the above information represents typical data for the product and should not be construed or used as an actual product performance specification.



For additional information and samples, call Ansell Healthcare at 800-952-9916 or visit us at [www.ansellhealthcare.com](http://www.ansellhealthcare.com)

Ansell Healthcare Products LLC | [ansellhealthcare.com](http://ansellhealthcare.com)  
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)  
Ansell, Derma Prene and IsoTouch are trademarks owned by Ansell Healthcare Products LLC or one of its affiliates. ©2009 All Rights Reserved.

SS330 04/2009

**Ansell**  
Everything you touch...we touch.