

Material 100% polychloroprene (neoprene); no natural rubber

Packaging 1 Sterile Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case

Product Numbers

SIZE (identified inside glove)	ASPN (Ansell Standard Product Number)
5 1/2	8511
6	8512
6 1/2	8513
7	8514
7 1/2	8515
8	8516
8 1/2	8517
9	8518

Dimensional Specifications (minimum length 295mm)

SIZE	Average Length mm (inches)	Average Palm Width mm (inches)
5 1/2	310 (12.2)	70 (2.8)
6	315 (12.4)	77 (3.0)
6 1/2	312 (12.3)	82 (3.2)
7	310 (12.2)	90 (3.5)
7 1/2	315 (12.4)	95 (3.7)
8	310 (12.2)	102 (4.0)
8 1/2	308 (12.1)	109 (4.3)
9	310 (12.2)	113 (4.4)

Thickness (single)

Data from sampling	Average Thickness
Fingertip, mm (mils)	0.20 (7.9)
Palm, mm (mils)	0.19 (7.5)
Cuff, mm (mils)	0.16 (6.3)

Physical Properties (ASTM Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level I)	1 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	22 (3,191)	28 (4,061)
Elongation (%) minimum	997	815
Modulus, MPa (p.s.i.) maximum	2.3 (334)	3.6 (522)
Force at Break (N)	13.1	15.5

Manufacturing specifications meet or exceed all the above.

Color	Green
Accelerators	None used in manufacturing process
Powder Free	Eliminates powder-related post-operative complications and sensitivity
Certification	Manufacturing Facility Registered to ISO 9001, EN46001 and ISO 13485 Tested for use with chemotherapy drugs (Contact Customer Technical Service at 888-HAND-ONE for information on test procedures and on the specific compounds tested.)
Quality Control	In compliance with ASTM D 3577-00, Type II
Sterilization	Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000
References	ANSI/ASQ Z1.4 Sampling Procedures and Tables for Inspection by Attributes
Country of Origin	Sri Lanka
Storage	Cool environment, shielded from ozone and UV light

*Acceptable Quality Level

For additional information and samples, call Ansell at 800-952-9916 or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
Ansell and Derma Prene are registered trademarks owned by Ansell Limited or one of its affiliates. ©2011 Ansell Healthcare Products LLC. All Rights Reserved.

