

Encore® Underglove

Encore® Underglove latex, powder-free surgical gloves offer exceptional dexterity, tactile sensitivity, outstanding strength and barrier protection when double-gloving is required. The Encore® Underglove is also designed to retain moisture and rehydrate your skin, despite continual glove-wearing and frequent contact with astringent hand cleaners and disinfectants. The Encore® Underglove can be used as an underglove, or for individual use.

Safety & Protection

Infection Control

Viral Penetration. Although not an FDA requirement, Ansell's premium line of surgical gloves is tested to ensure protection against viral permeation. Encore® has been tested and passed the viral penetration according to ASTM F1671-97a.

Cuff Roll-Down. This glove has Ansell's patent-pending SureFit Technology™, which provides surgeons, nurses and other healthcare professionals with superior barrier protection during operating room procedures.

Bone Cement. As long as the preparation and use directions on the bone cement are followed, this glove provides a protective barrier for handling mixed bone-cement product.

Embryo Use. This glove is appropriate and accepted for embryo use.

Skin Protection

Low Protein Content. This glove contains 50 µg/g or less of water/extractable protein per gram.

Caution: Safe use of these gloves by latex-sensitized individuals has not been established.

Moisturizing. Encore® Underglove features Hydrasoft Technology® which utilizes two moisturizers: skin-friendly glycerine, which attracts and retains water, and dimethicone, which prevents the skin from drying out.



Comfort & Fit

Donning

Alcohol Resistant. Ansell's proprietary polyurethane inner coating will not break down. It is resistant to alcohol-based hand disinfectants and adds another layer of safety between your skin and the latex.

Damp-Hand Donnability. Polymer coating for dry and damp-hand donnability.

Double Donnability. Thickness allows for easy and comfortable double-gloving. Can be used as an underglove, or for individual use.

Hand Fatigue

Shape of Glove/Former. Designed with Anatomic Curved Former for a more natural fit and less fatigue. This exclusive Ansell former is the result of a staff dedicated to the development and continuous improvement of formers.

Stretchability (Modulus). Because this glove is latex, wearers experience a more secure fit that provides greater comfort without being baggy.

Tactile Sensitivity

Dependability. Micro-thin design provides soft comfort and allows for fingertip sensitivity during double gloving or individual use.

Grip

Application Driven. Anti-slip smooth finish prevents slip.

ENCORE®
UNDERGLOVE



Ansell

Everything you touch...we touch.