

**Material** Natural rubber from latex in compliance with ASTM D 1076

**Packaging** 1 Sterile Pair/Polyethylene Package – 50 Pairs/Dispenser Box – 4 Boxes/Case

### Product Numbers

SIZE (identified inside glove)	ASPN (Ansell Standard Product Number)
5 1/2	2018455
6	2018460
6 1/2	2018465
7	2018470
7 1/2	2018475
8	2018480
8 1/2	2018485
9	2018490

### Dimensional Specifications (minimum length 295mm)

SIZE	Average Length mm (inches), ±4mm	Average Palm Width mm (inches)
5 1/2	296 (11.7)	70 (2.7)
6	297 (11.7)	76 (3.0)
6 1/2	298 (11.7)	82 (3.2)
7	298 (11.7)	89 (3.5)
7 1/2	298 (11.7)	95 (3.7)
8	300 (11.8)	102 (4.0)
8 1/2	300 (11.8)	108 (4.3)
9	303 (11.9)	115 (4.5)

### Thickness (single)

Data from sampling	Average Thickness	Range
Fingertip, mm (mils)	0.18 (7.1)	±0.018 (0.7)
Cuff, mm (mils)	0.11 (5.0)	±0.011 (0.5)

### Physical Properties (ASTM, EN and ISO Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level GI)	1 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	24 (3,480)	18 (2,610)
Elongation (%) minimum	750	560
Modulus, MPa (p.s.i.) maximum	5.5 (800)	N/A

*Manufacturing specifications meet or exceed all the above.*

<b>Color</b>	Teal
<b>Accelerators</b>	Zinc dibutyldithiocarbamate (ZDBC) Zinc diethyldithiocarbamate (ZDEC)
<b>Surface Treatment</b>	External: chlorinated and siliconized Internal: polyurethane donning coating with moisturizer
<b>Powder Free</b>	Eliminates powder-related post-operative complications and sensitivity
<b>Certification</b>	Manufacturing Facility Registered to ISO 9001, 13485 and CE certification Contains 50 micrograms or less of total water-extractable protein per gram
<b>Quality Control</b>	In compliance with ASTM D 3577-06, Type I
<b>Sterilization</b>	Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000
<b>References</b>	ANSI/ASQ Z1.4 Sampling Procedures and Tables for Inspection by Attributes
<b>Country of Origin</b>	Sri Lanka
<b>Storage</b>	Cool environment, shielded from ozone and UV light

\*Acceptable Quality Level

*Except for applicable ASTM requirements, the above information represents typical data for the product and should not be construed or used as an actual product performance specification.*

**For additional information and samples, call Ansell at 800-952-9916 or visit us at [www.ansellhealthcare.com](http://www.ansellhealthcare.com)**

Ansell Healthcare Products LLC | [ansellhealthcare.com](http://ansellhealthcare.com)  
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)

Ansell and Encore are registered trademarks owned by Ansell Healthcare Products LLC or one of its affiliates.  
©2011 Ansell Healthcare Products LLC. All Rights Reserved.

SS310 11/2011



**Ansell**  
Everything you touch...we touch.