

Technical Data

No Powder™ Latex, Powder-Free Surgical Gloves, 860x

Material Natural rubber latex

Packaging 1 Sterile Pair/Package – 50 Pairs/Dispenser Box – 4 Boxes/Case (200 Pairs/Case)

Product Numbers

SIZE (identified inside glove)	ASPN (Ansell Standard Product Number)
5 1/2	N/A
6	8602
6 1/2	8603
7	8604
7 1/2	8605
8	8606
8 1/2	8607
9	8608

Dimensional Specifications (minimum length 295mm)

SIZE	Average Length mm (inches)	Average Palm Width mm (inches)
5 1/2	N/A	N/A
6	300 (11.8)	80 (3.15)
6 1/2	300 (11.8)	85 (3.35)
7	300 (11.8)	90 (3.54)
7 1/2	300 (11.8)	95 (3.74)
8	300 (11.8)	100 (3.94)
8 1/2	300 (11.8)	110 (4.33)
9	300 (11.8)	115 (4.53)

Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.23 (9.1)
Palm, mm (mils)	0.22 (8.7)
Cuff, mm (mils)	0.25 (9.8)

Physical Properties (ASTM, EN, and ISO Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level I)	1 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	24 (3,480)	18 (2,610)
Elongation (%) minimum	750	560
Modulus, MPa (p.s.i.) maximum	5.5 (800)	N/A
Force at Break (N)	20	20

Ansell manufacturing specifications meet or exceed all the above.

Color	Cream
Accelerators	Zinc diethyldithiocarbamate (ZDEC) Zinc 2-mercaptobenzothiazole
Powder Free	Eliminates powder-related post-operative complications and sensitivity
Donning Agent	None (powder-free)
Residual Powder Content	<2 mg/glove (typical range, 0.10–0.20 mg/glove)
Certification	Manufacturing facility registered to ISO 9001 and ISO 13485
Quality Control	In compliance with ASTM D 3577, Type 1
Sterilization	Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000
References	ANSI/ASQC Z1.4 Sampling Procedures and Tables for Inspection by Attributes
Country of Origin	Malaysia
Storage	Cool environment, shielded from ozone and UV light

*Acceptable Quality Level

Except for applicable ASTM requirements, the above information represents typical data for the product and should not be construed or used as an actual product performance specification.



For additional information and samples, call Ansell Healthcare at 800-952-9916 or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC | ansellhealthcare.com
200 Schulz Drive | Red Bank, NJ 07701 U.S.A. | Phone: (800) 952-9916 (U.S.A. only)
Ansell and No Powder are trademarks owned by Ansell Limited or one of its affiliates. ©2009 All Rights Reserved.

SS340 04/2009

Ansell
Everything you touch...we touch.